


# NeuroSafe


How can technology of local brain cooling assist World Healthcare





# What is craniocerebral hypothermia?

Craniocerebral hypothermia is controlled cooling of the brain, which reduces metabolism of neurons, enhances strength of cell membranes, decreases their requirement in oxygen, and as a result makes them more viable








# When is craniocerebral hypothermia necessary?


In case of diseases which are accompanied by overheating of the brain and result in destruction of the neurons, causing death and disability

# Wide range of application






Why is craniocerebral hypothermia necessary?




# Priorities of World Healthcare





# The solution



# NeuroSafe

Compact, portable model: beta testing, registration process is progress

Compressor

Cooling of one patient

Weight (with coolant) 12 кг

Cooling of the brain T0 no lower than 29°C

Speed of cooling - 1 °C за 45 минут

Time of cooling of coolant from +25°C to 0°C: 20 min.

Display with user-friendly interface

Easy transportation

Automatic protocol of the procedure

Engineering - ENCATA

Prototype – June 2021

Registration – September 2021





# THE-01

## MVP (Original model)

Registered: PY ФСР 2011/11788 от 12.12.2017, Declaration of compliance РОССТ РУ Д-РУ.ИМ25.В.00024/18 (30.10.2018)

Serial Production phase: «Izhevsk Mechanical Factory» Concern «Kalashnikov»

10 Patents

Cooling of the brain  $T_0$  no lower than  $29^{\circ}\text{C}$

Cooling of body  $T^0$   $36,6 \dots 33^{\circ}\text{C}$

Time of cooling of coolant from to  $0^{\circ}\text{C}$  - 30 min

Speed of cooling -  $1^{\circ}\text{C}$  per 45 minutes

Cooling of two patients simultaneously

Automatic protocol of the procedure

3 channels of monitoring of patients body  $T^0$

Display with user-friendly interface

Automatic procedure control

Size  $94 \times 45 \times 52$  cm


Weight (with coolant) 70 kg

Reusable cryoapplicators





# The effect



Improves pharmacoeconomical characteristics of the medical care

Improves cognitive abilities and reduces fatigue

Reduces time in the resuscitation units by 1,5 times

Halves neurological deficit

Halves death rate among patients in critical conditions

More than **5 000** people

have been treated in Russia, Kazakhstan and Uzbekistan clinics with defined therapeutic effect and without side effects and complications






# How does it work?


«When craniocerebral hypothermia initiates neuroprotection,  
one can benefit from it long after cooling»\*

*\*Nina M. Rzechorzek et al (2016) Hypothermic Preconditioning Reverses Tau Ontogenesis in Human Cortical Neurons and is Mimicked by Protein Phosphatase 2A Inhibition*



In the area of lesion **E** edema is formed, which in combination with central heat inflows **D** does not allow cold blood flow to reach core of hyperthermia area, Temperature rises up to 41°C, which leads to damage and destruction of neurons

Brain T°C in case of stroke before using «THE-01»



Brain T°C after using «THE-01» after 4 hours of the procedure

# Software

Allows to control the procedure

Allows to monitor temperature of the patient and coolant during the procedure

Tc - coolant T°C

Th - helmet (cryoapplicator) T°C

Tt - tympanic T°C (measured in the ear area)

Ta - axillary T°C (measured in the armpit)

Allows to protocol the procedure

Allows to adjust settings for different nosologies

**Preparing equipment for the procedure**  
When Tc (coolant T°) reaches +5°C,  
THE is ready for the procedure

Set Tc **Tc +5°C**

Th: +23,8°C  
Tt: +23,9°C  
Ta: +27,1°C

Settings  
Length of the procedure: 00:01  
Th: 05,0°C  
Tt: 15,0°C    Ta: 15,0°C

**Start**

Protocol    Settings

# Software Settings

Allows to set temperature of the coolant

Allows to set length of the procedure

Allows to set goal temperature of the patient and coolant during the procedure

Th – helmet (cryoapplicator) T°C

Tt – tympanic T°C (measured in the ear area)

Ta – axillary T°C (measured in the armpit)

THE is ready to cool

Set Tc **0 °C**

Th: +23,9°C

Tt: +24 °C

Ta: +24 °C

Settings

Length of the procedure: **48:00**

Th: **13 : 5**

Tt: **33:0**      Ta: **33:0**

**Start**

**Protocol**      **Settings**

The screenshot shows a software interface for temperature control. At the top, a green box displays 'THE is ready to cool'. Below this, the coolant temperature is set to 'Set Tc 0 °C'. A central panel shows a human figure with three temperature measurement points: 'Th: +23,9°C' (helmet), 'Tt: +24 °C' (tympanic), and 'Ta: +24 °C' (axillary). A 'Settings' section at the bottom allows for adjusting the 'Length of the procedure' to '48:00', and setting target temperatures for 'Th: 13 : 5', 'Tt: 33:0', and 'Ta: 33:0'. A 'Start' button with a green progress bar is visible, along with 'Protocol' and 'Settings' buttons.

# Protocol of the procedure

Patients medical details

Details of CHH procedure and patient's parameters before the treatment

Length of the Procedure

**Preparing equipment for the procedure**

When Tc (coolant T°) reaches + 5°C,  
THE is ready for the procedure

Name of Patient №1

Date: \_\_\_\_\_  
Gender: \_\_\_\_\_ Age: \_\_\_\_\_  
Diagnosis: \_\_\_\_\_

Criteria of the inclusion:  
Start of CCH: \_\_\_\_\_  
Finish of CCH: \_\_\_\_\_  
Artificial ventilation: \_\_\_\_\_  
Sedation/therapy: \_\_\_\_\_  
NIHSS: \_\_\_\_\_  
Glasgo: \_\_\_\_\_  
Body T: \_\_\_\_\_ Tympanic (ear) T: \_\_\_\_\_  
Blood pressure: \_\_\_\_\_ Heart rate: \_\_\_\_\_


0,5 hour  
1 hour  
2 hours  
4 hours  
8 hours  
16 hours  
20 hours  
24 hours



# Data Analysis

Visual representation of temperature of the patient and helmet during the procedure

- Th – helmet (cryoapplicator) T°C
- Tt – tympanic T°C (measured in the ear area)
- Ta – axillary T°C (measured in the armpit)





|   | Document              | Date of publication | Expiry date |
|---|-----------------------|---------------------|-------------|
| Device for human body local cooling   | RF Patent no. 83369   | 16.10.2012          | 30.03.2026  |
| Device for therapeutic hypothermia induction                                    | RF Patent no. 126262  | 27.03.2013          | 28.08.2022  |
| Device for correction of cerebral hyperthermia                                  | RF Patent no. 2615283 | 04.04.2017          | 12.11.2035  |
| «Method of diagnostics and management of cerebral hypothermia syndrome»         | RF Patent no. 2645927 | 15.06.2017          | 15.06.2037  |
| «Device for conduction of controlled therapeutic hypothermia of the brain»      | RF Patent no. 2653794 | 14.05.2018          | 25.08.2037  |
| Application for «Device for induction of controlled hypothermia of the brain»   | PCT/RU2017/000633     | In progress         | In progress |
| Application for «Method of therapy of knee joints using deep local hypothermia» | RF Patent no. 2653794 | 14.05.2018          | 25.08.2037  |

# Why us?

Individual approach and widest range of application:  
suitable for conscious state, coma, medical anesthesia

Safest and most convenient method:  
no pharmacological assistance required, no complications, no side effects

Reduction of expenses to hospitals

The most effective hypothermia device on the market,  
which has the widest range of application

The most attractive price for equipment and consumables on the market  
cryoapplicators – reusable consumables (6 months)

**N°1**  
**in the world**  
Price/Quality ratio parameter



# Reasons to consider us

## Healthcare problems in EU

- Stroke is the leading cause of death in EU
- Europe has highest proportion of heavy episodic drinking in the World
- Over 66 million people have chronic depression (7% of population)
- Almost 74 million people suffer from migraines (10% of population)

Market volume in USA, Canada, EU, Japan – 250 000 units

Export potential: annual market growth in developed countries (USA, Canada, EU, Japan) 7-10%

# €5 Billion EUR

Market volume in USA, Canada, EU, Japan

## 37%



Of French ICUs are equipped with general hypothermia equipment



## 63%

Of French ICUs are not equipped with general hypothermia equipment

# NeuroSafe



## What we do:

Development and production of craniocerebral therapeutic hypothermia equipment line

Introduction of craniocerebral hypothermia into medical practice

Development of methodologies of craniocerebral hypothermia

Scientific research of craniocerebral hypothermia



## Mission statement:

Development of effective and economically efficient hypothermia technologies with aim to improve resuscitation and rehabilitation outcomes:

- reduction of mortality
- reduction of rehabilitation period
- protection of health and preventive healthcare
- improvement of quality of life



## Support:

Agency of Strategic Initiatives (RUS)

Skolkovo Innovation Centre (RUS)

# We are ready for cooperation!

Galina Repina

Leading Marketing Specialist

Mob. +7 915 421 13 50

E-mail: [neurosafe@mail.ru](mailto:neurosafe@mail.ru)

Сайт: [www.neurosafe.ru](http://www.neurosafe.ru)

